

Human IgG1 Isotype Control (5Vk2/37.7H) mAb

Catalog No.	FHJ92810A FHJ92810B	Quantity:	100 µg 1.0 mg
Description:	Human monoclonal IgG1 kappa antibody, purified by affinity chromatography, with no known specificity, for use as an isotype control.		
Target:	No known specificity		
Concentration:	~ 1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Clone:	5Vk2/37.7H		
Purity:	> 95% determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	FCM, Isotype control, In vivo control		
Storage & Stability:	Store at 2-8°C for immediate use (1-2 weeks) or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

