

SLITRK6

Humanized Anti-Human SLITRK6 (SAA0827) mAb

Catalog No.	FHK06810A	Quantity:	50 µg
	FHK06810B		100 µg
	FHK06810C		1.0 mg

Alternate Names: SLIT and NTRK-like protein 6

Description: SLITRK6 is member of the SLITRK family of integral membrane proteins with 2 N-terminal leucine-rich repeat (LRR) domains. Most SLITRK members also have C-terminal regions that share homology with neutrotrophin receptors and are expressed predominately in neurite tissues and have neurite-modulating activity. An antibody-drug conjugate (ADC) is in clinical trials for urothelial cancer.

UniProt ID (target): Q9H5Y7

Target: Human SLITRK6

Concentration: 1.0 mg/ml, lot specific

Source: XtenCHO

Isotype: Human IgG2 kappa

Formulation: 0.01M PBS buffer, pH 7.4

Purity: > 95% by SDS-PAGE

Purification: Protein A/G affinity chromatography

Applications: ELISA, FCM, IF

Storage & Stability: Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year.
Avoid freeze/thaw cycles.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.