


Human IgG1 (LALA P329G) Isotype Control (13R4) mAb

Catalog No.	FHK17820A FHK17820B	Quantity:	100 µg 100 µg
Description:	Human IgG1 (LALA P329G) monoclonal antibody (13R4) for use as an <i>in vivo</i> control or isotype control.		
Target:	No known specificity		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG1 (LALA P329G)		
Formulation:	PBS, pH 7.4 containing 0.2 M Arg		
Purity:	> 95% by SDS-PAGE		
Clone:	13R4		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	FCM, Isotype Control, In vivo Control		
Manufacturer:	Antibody System		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

