

TNFRSF25

Mouse Anti-Human TNFRSF25 (PTX25) mAb

Catalog No.	FHK23210A FHK23210B	Quantity:	50 µg 100 µg
Alternate Names:	Tumor necrosis factor receptor superfamily member 25, Apo-3, Apoptosis-inducing receptor AIR, Apoptosis-mediating receptor DR3, Apoptosis-mediating receptor TRAMP, Death receptor 3		
Description:	TNFRSF25 is the receptor for TNFSF12/APO3L/TWEAK. This receptor is expressed preferentially in the tissues enriched in lymphocytes, and it may play a role in regulating lymphocyte homeostasis. This receptor has been shown to stimulate NF-kappa B activity and regulate cell apoptosis. The signal transduction of this receptor is mediated by various death domain containing adaptor proteins.		
UniProt ID (target):	Q93038		
Target:	Human TNFRSF25		
Concentration:	1.0 mg/ml, lot specific		
Clone:	PTX25		
Isotype:	Mouse IgG2 lambda		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

