

GPC2

Mouse Anti-Human GPC2/Glypican-2 (D3) mAb

Catalog No.	FHK80620A FHK80620B	Quantity:	50 µg 100 µg
Alternate Names:	Glypican-2, Cerebroglycan		
Description:	Glypican-2 (cerebroglycan) is a cell surface proteoglycan bearing heparan sulfate, that participates in cell adhesion and is thought to regulate the growth and guidance of axons. GPC2 is highly expressed in neuroblastoma and is considered a potential therapeutic target for this disease.		
UniProt ID (target):	Q8N158		
Target:	Human Glypican-2		
Concentration:	1.0 mg/ml, lot specific		
Clone:	D3		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE and SEC-HPLC		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

