

BTN2A1

Mouse Anti-Human BTN2A1 (SAA2136) mAb, APC

Catalog No.	FHN38913A FHN38913B	Quantity:	50 tests 100 tests
Alternate Names:	Butyrophilin subfamily 2 member A1, BTN2A1		
Description:	Mouse Anti-human BTN2A1 (SAA2136) was affinity purified from cell culture supernatant, targets Human BTN2A1, followed by conjugation to allophycocyanin.		
UniProt ID (target):	Q7KYR7		
Target:	Human Butyrophilin subfamily 2 member A1		
Concentration:	0.5 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a		
Formulation:	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.		
Purity:	> 95% by SDS-PAGE		
Clone:	SAA2136		
Purification:	Protein A/G affinity chromatography , ion exchange chromatography		
Applications:	FCM, ELISA		
Conjugate:	Allophycocyanin (APC)		
Storage & Stability:	Store at 2-8°C for up to 1 year. Protect from light. Do not freeze		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

