

Mouse IgG2a, kappa Isotype Control (MOPC-173) mAb, FITC

Catalog No.	FMB96751A	Quantity:	50 tests
	FMB96751B		100 tests

Description: Mouse monoclonal IgG2a, kappa antibody, purified by affinity chromatography, conjugated to fluorescein isothiocyanate (FITC).
Useful as an isotype control, with no known specificity.

Target: No known specificity

Concentration: 0.5 mg/ml, lot specific

Source: Mammalian cell line

Isotype: Mouse IgG2a, kappa

Formulation: 0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300

Clone: MOPC-173

Purity: > 95% determined by SDS-PAGE

Purification: Protein A/G affinity chromatography

Conjugate: FITC

Applications: FCM, Isotype control, In vivo control

Storage & Stability: Store at 2-8°C for up to 1 year. Protect from light. **Do not freeze.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.