

Kel

Rat Anti-Mouse CD238/KEL (MIMA-8) mAb

Catalog No.	FMD56310A FMD56310B	Quantity:	100 µg 1.0 mg
Alternate Names:	Kell blood group glycoprotein homolog, CD238		
Description:	Affinity purified Rat Anti-Mouse CD238/KEL (MIMA-8) targets mouse Kell blood group glycoprotein homolog		
UniProt ID (target):	Q9EQF2		
Specificity:	Mouse Kell blood group glycoprotein homolog Cross reacts with human CD238/KEL		
Source:	Mammalian cell line		
Isotype:	Rat IgG1 kappa		
Clone:	MIMA-8		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	Flow cytometry		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

