

Adam17

Rat Anti-Mouse CD156b/ADAM17 (A12) mAb

Catalog No.	FMF68810A FMF68810B	Quantity:	100 µg 1.0 mg
Alternate Names:	Disintegrin and metalloproteinase domain-containing protein 17, ADAM 17, CD156b, TNF-alpha convertase, TNF-alpha-converting enzyme		
Description:	Rat Anti-Mouse CD156b/ADAM17 (A12) targets human and mouse TNF-alpha-converting enzyme (ADAM 17, CD156b)		
UniProt ID (target):	Q9Z0F8		
Specificity:	Human, Mouse TNF-alpha-converting enzyme (ADAM 17, CD156b)		
Source:	Mammalian cell line		
Isotype:	Rat IgG1 kappa		
Clone:	A12		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	ELISA, Flow cytometry		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.