

Cd276

Rat Anti-Mouse/Human CD276 / B7-H3 (m276) mAb

Catalog No.	FMJ04020A FMJ04020B	Quantity:	50 µg 100 µg
Alternate Names:	CD276 antigen, B7 homolog 3, B7-H3, Costimulatory molecule, CD276		
Description:	Clone m276 is a Rat IgG1 monoclonal antibody that is cross-species reactive with mouse, human and rat CD276.		
UniProt ID (target):	Q8VE98		
Target:	Mouse/Human (Rat) CD276 / B7-H3		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Rat IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

