

Cd248

Human Anti-Mouse CD248 (G78) mAb, PE

Catalog No.	FMJ61012A FMJ61012B	Quantity:	50 tests 100 tests
Alternate Names:	Endosialin, Tumor endothelial marker 1, CD248		
Description:	Clone G78 is a Human IgG1 lambda monoclonal antibody that is cross-species reactive with mouse, human CD248, purified by affinity chromatography and conjugated to phycoerythrin.		
UniProt ID (target):	Q91V98		
Target:	Human, Mouse CD248		
Concentration:	0.5 mg/ml, lot specific		
Conjugate:	Phycoerythrin (PE)		
Isotype:	Human IgG1 lambda		
Formulation:	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Clone:	G78		
Purification:	Protein A/G affinity chromatography		
Applications:	ELISA, FCM (0.5µg/test)		
Storage & Stability:	Store at 2-8°C for 1 year. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com