

Mouse IgG1, kappa Isotype Control (MOPC-21) mAb

Catalog No.	FMJ92820A FMJ92820B	Quantity:	50 µg 100 µg
Description:	Mouse monoclonal IgG1, kappa antibody, purified by affinity chromatography, with no known specificity, for use as an isotype control.		
Target:	No known specificity		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Mouse IgG1, kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Clone:	MOPC-21		
Purity:	> 95% determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM, Isotype control		
Storage & Stability:	Store at 2-8°C for immediate use (1-2 weeks) or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.