

Mouse IgG3 Isotype Control (B1-8) mAb

Catalog No.	FMJ93040A FMJ93040B	Quantity:	100 µg 1.0 mg
Description:	Mouse IgG3 monoclonal antibody is useful as an isotype control.		
Specificity:	No known specificity		
Isotype:	Mouse IgG3		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM, Control		
Storage & Stability:	Stable at 2-8°C for 1 week or for up to 1 year at -20°C to -80°C. It is recommended to prepare working aliquots of undiluted product and store -20°C to -80°C. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

