

DII3

Mouse Anti-Human DLL3 (SC16.4) mAb

Catalog No.	FMK06410A FMK06410B	Quantity:	50 µg 100 µg
Alternate Names:	Delta-like protein 3, Drosophila Delta homolog 3, Delta3, M-Delta-3		
Description:	Mouse Anti-Human DLL3 (SC16.4) mAb reacts with Human, Mouse, Rat Delta-like protein 3.		
UniProt ID:	O88516		
Specificity:	Recognizes Human, Mouse, Rat Delta-like protein 3		
Isotype:	Mouse IgG2a kappa		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Purification:	Protein A/G affinity chromatography		
Applications:	ELISA, FCM		
Storage & Stability:	Stable at 2-8°C for 1 week or for up to 1 year at -20°C to -80°C. It is recommended to prepare working aliquots of undiluted product and store -20°C to -80°C. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

