

Canine IgG4, kappa Isotype Control (HyHEL-10) mAb

Catalog No.	FQK18110A FQK18110B	Quantity:	50 µg 100 µg
Description:	Canine IgG4 kappa monoclonal antibody is useful as an isotype control in flow cytometry.		
Specificity:	No known specificity		
Isotype:	Canine IgG4, kappa light chain		
Clonality:	Monoclonal		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4 containing 0.09% sodium azide		
Purification:	Protein A/G affinity chromatography		
Purity:	> 95% by SDS-PAGE		
Applications:	FCM, isotype control		
Storage & Stability:	Stable at 2-8°C for 1 week or for up to 1 year at -20°C to -80°C. It is recommended to store aliquots of the stock solution at -20°C to -80°C. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.