

Canine IgG4, lambda Isotype Control (HyHEL-10) mAb

Catalog No.	FQK18120A FQK18120B	Quantity:	50 µg 100 µg
Description:	Recombinant Canine IgG4 lambda monoclonal antibody with no known specificity, purified by Protein A/G affinity chromatography, for use as an isotype control.		
UniProt ID:	P00698		
Specificity:	No known specificity		
Source:	Mammalian cell line		
Isotype:	Canine IgG4 lambda		
Clone:	HyHEL-10		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4, 0.09% sodium azide		
Purification:	Protein A or G affinity chromatography		
Purity:	>95% as determined by SDS-PAGE		
Applications:	FCM, Isotype control		
Storage & Stability:	Stable at 2-8°C for 1 week or for up to 1 year at -20°C to -80°C. It is recommended to prepare working aliquots of undiluted product and store -20°C to -80°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.