


THPO

Anti-Romiplostim (1A430) mAb

Catalog No.	MHE26001A MHE26001B	Quantity:	50 µg 100 µg
Alternate Names:	Thrombopoietin, Romiplostim, Nplate, Romiplate, Megakaryocyte colony-stimulating factor, Megakaryocyte growth and development factor, MGDF, Myeloproliferative leukemia virus oncogene ligand, THPO		
Description:	Affinity purified mouse Anti-Romiplostim monoclonal antibody (clone 1A430) is useful as the capture antibody for Romiplostim ELISA.		
UniProt ID (target):	P40225		
Specificity:	Romiplostim (human)		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Mouse IgG2a		
Formulation:	Sterile filtered PBS, 0.2M Arg, pH 6.0		
Immunogen:	Recombinant Romiplostim protein		
Purity:	> 95% by SDS-PAGE and SEC-HPLC		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Clone:	1A430		
Applications:	ELISA capture antibody, use with detection antibody MHE26002.		
Manufacturer:	Antibody System		
Storage & Stability:	Store at 2-8°C for immediate use (1 week). Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

