

THPO

Anti-Romiplostim (1A432) mAb

Catalog No.	MHE26003A MHE26003B	Quantity:	50 µg 100 µg
Alternate Names:	Thrombopoietin, Romiplostim, Nplate, Romiplate, Megakaryocyte colony-stimulating factor, Megakaryocyte growth and development factor, MGDF, Myeloproliferative leukemia virus oncogene ligand, THPO		
Description:	Affinity purified mouse Anti-Romiplostim monoclonal antibody (clone 1A432) was expressed in mammalian cells.		
UniProt ID (target):	P40225		
Specificity:	Romiplostim (human)		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Mouse IgG1		
Formulation:	Sterile filtered 0.01M PBS, pH 7.4		
Immunogen:	Recombinant Romiplostim protein		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Clone:	1A432		
Applications:	ELISA		
Storage & Stability:	Store at 2-8°C for immediate use (1 week). Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com