

Nrn1

Mouse Anti-Murine NRN1/Neuritin (1A440) mAb

Catalog No.	MMN45401A MMN45401B	Quantity:	50 µg 100 µg
Alternate Names:	Neuritin, Candidate plasticity gene 15 protein		
Description:	Recombinant antibody mouse anti-murine Neuritin (1A440) was expressed in a mammalian cell line and purified using Protein A/G affinity chromatography.		
UniProt ID (target):	Q8CFV4		
Target:	Mouse Neuritin		
Concentration:	~1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Mouse IgG2b		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	1A440		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	ELISA, Western blot		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com