

CD8A

Mouse Anti-Human CD8a (Clone FK18) mAb

Catalog No.	MON1008-1 MON1008-5	Quantity:	1 ml 5 ml
Alternate Names:	T-cell surface glycoprotein CD8 alpha chain, T-lymphocyte differentiation antigen T8/Leu-2		
Description:	MON1008 reacts with CD8, the 33 kDa T8 antigen present on human cytotoxic/suppressor T-cell subset. It strongly blocks the cytolytic activity of T8-positive CTLs. It also reacts with a subpopulation of Rhesus monkey lymphocytes which is identical to the population of human T8 positive cells.		
UniProt ID:	P01732		
Specificity:	CD8		
Isotype:	Mouse IgG3		
Clone:	FK18		
Formulation:	PBS with 0.02% sodium azide		
Applications:	IF: suggested dilution 1:10 in PBS, 1% BSA, 0.02% sodium azide		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

