

EGFR

Mouse Anti-Human EGFR/ErbB-1 (Clone 2E9) mAb

Catalog No.	MON1033-1 MON1033-5	Quantity:	1 ml 5 ml
Alternate Names:	Epidermal growth factor receptor, proto-oncogene c-ErbB-1, Receptor tyrosine-protein kinase erbB-1		
Description:	MON1033 reacts with a protein determinant of the extra cellular domain of the human EGFR, and does not cross react with EGFR from murine cells. The antibody shows binding competition with EGF. The antibody reacts in immunoprecipitation with the functional EFR protein-tyrosine kinase complex.		
Concentration:	0.1 mg/ml		
UniProt ID:	P00533		
Specificity:	Human Epidermal Growth Factor Receptor		
Host:	Mouse		
Isotype:	IgG1		
Clone:	2E9		
Formulation:	PBS containing 1% BSA and 0.1% sodium azide		
Applications:	IF, IP		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

