

## ITGA5

### Mouse Anti-Human Fibronectin Receptor / VLA-5 (Clone NKI-SAM1) mAb

<b>Catalog No.</b>	MON1034	<b>Quantity:</b>	1 ml
<b>Alternate Names:</b>	Integrin alpha-5, CD49 antigen-like family member E, Fibronectin receptor subunit alpha, Integrin alpha-F, VLA-5, CD49e		
<b>Description:</b>	Integrin alpha-5/beta-1 (ITGA5:ITGB1) is a receptor for fibronectin and fibrinogen. It recognizes the sequence R-G-D in its ligands.		
<b>Concentration:</b>	100 µg/ml		
<b>UniProt ID:</b>	P08648		
<b>Specificity:</b>	Human Fibronectin receptor		
<b>Host:</b>	Mouse		
<b>Isotype:</b>	IgG2b		
<b>Clone:</b>	NKI-SAM1		
<b>Formulation:</b>	PBS containing 0.1% sodium azide		
<b>Applications:</b>	IHC-F		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

