

Mouse Anti-Human HLA-DR + DP (Clone HL-40) mAb

Catalog No.	MON1083	Quantity:	100 µg
Description:	MON1083 reacts with common epitope on beta-chain of human HLA-DR and HLA-DP. DR and DP are the isotypes of human MHC Class II molecules expressed on antigen-presenting cells (APC; dendritic cells, B lymphocytes, monocytes, macrophages)		
Concentration:	1 mg/ml		
Specificity:	Human and porcine HLA-DR + DP		
Host:	Mouse		
Immunogen:	Raji Burkitt's lymphoma cell line		
Isotype:	IgG2a		
Formulation:	PBS, pH 7.4 containing 15 mM sodium azide		
Cross-Reactivity:	Pig		
Applications:	Flow cytometry: suggested concentration 1 µg/ml Immunoprecipitation IHC-P		
Storage & Stability:	Store antibody at 2-8°C until expiry date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

