

CD44

Mouse Anti-Human CD44 (Clone MEM-85) mAb, FITC

Catalog No.	MON1123F	Quantity:	0.50 ml
Alternate Names:	CD44 antigen, CDw44, Epican, Extracellular matrix receptor III, ECMR-III, GP90 lymphocyte homing/adhesion receptor, HUTCH-I, Heparan sulfate proteoglycan, Hermes antigen, Hyaluronate receptor, Phagocytic glycoprotein I, PGP-I, CD44		
Description:	MON 1123 reacts with both cell surface-expressed and soluble form of CD44 antigen (Phagocyte glycoprotein 1), a 80-95 kDa transmembrane glycoprotein (hyaladherin family) present on the most of cells and tissues (leukocytes, endothelial cells, mesenchymal cells, etc.); it is negative on platelets and hepatocytes.		
Concentration:	1 mg/ml		
UniProt ID:	P16070		
Specificity:	CD44		
Host:	Mouse		
Immunogen:	Leukocytes of a patient with LGL Type Leukemia.		
Isotype:	IgG2b		
Clone:	MEM-85		
Conjugate:	FITC		
Formulation:	PBS, pH 7.4 containing 4% BSA and 0.1% sodium azide		
Applications:	FC		
Storage & Stability:	Store undiluted antibody at 2-8°C until expiration date. Avoid direct light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

