

CD99

Mouse Anti-human CD99 (Clone MEM-131) mAb

Catalog No.	MON1158	Quantity:	100 µg
Alternate Names:	CD99 antigen, 12E7, E2 antigen, Protein MIC2, T-cell surface glycoprotein E2		
Description:	MON 1158 reacts with CD99R, an epitope restricted to a subset of CD99 molecule expressed on myeloid cells, NK cells and T lymphocytes.		
Concentration:	1 mg/ml		
UniProt ID:	P14209		
Specificity:	CD99		
Host:	Mouse		
Immunogen:	HPB-ALL human peripheral blood leukemia T-cell line		
Isotype:	IgM		
Clone:	MEM-131		
Formulation:	PBS containing 15 mM sodium azide		
Applications:	WB, FC, IP, IHC-P		
Storage & Stability:	Store undiluted antibody at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

