

CGB3

Mouse Anti-Human Beta-HCG (Clone 2092) mAb

Catalog No.	MON2001-1 MON2001-5	Quantity:	1.0 ml 5.0 ml
Alternate Names:	Choriogonadotropin subunit beta 3, CG-beta, Beta-hCG		
Description:	MON 2001 reacts with high specificity with human chorionic gonadotrophin beta		
Concentration:	20 µg/ml		
UniProt ID:	P0DN86		
Specificity:	beta-HCG		
Host:	Mouse		
Immunogen:	Human chorionic gonadotropin		
Isotype:	IgG1		
Clone:	2092		
Formulation:	Tissue culture supernatant with PBS, 1% BSA, 10mM sodium azide		
Cross-Reactivity:	FSH: not detectable MSH: not detectable LH: slightly detectable		
Applications:	IHC: suggested dilution 1:5 - 1:25 IF		
Storage & Stability:	Store at 2-8°C for up to 1 year		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

