

## Staphylococcus spp, clone 17-5 Monoclonal Antibody

**Catalog No.:** MON2008-1

**Quantity:** 1 ml

**Specificity**

The antibody reacts with several *Staphylococcus* strains using ELISA. It does not react with several members of the *Enterobacteriaceae*.

**Immunoglobulin type**

Murine IgG<sub>3</sub>

**Use**

The antibody is useful for immunofluorescence, ELISA and agglutination test.

**Instructions for use**

For use in immunofluorescence the antibody should be used undiluted, for ELISA dilute 10x. For use in agglutination test prepare bacteria suspension in a physiological salt solution and place 10 µl on glass slide. Add 10 µl antibody solution, mix well and watch clumping. Incorporate always negative and positive controls.

**Presentation**

1 ml diluted ascites in PBS with approximately 23 µg/ml antibody and 0.1% Na-azide and 1% BSA, sufficient for 100-200 tests. Cat. nr. MON2008-1. The product is sterile filtered.  
5 ml diluted ascites in PBS with approximately 23 µg/ml antibody and 0.1% Na-azide and 1% BSA, sufficient for 500-1000 tests. Cat. nr. MON2008-5. The product is sterile filtered.

**FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.**



**Cell Sciences, Inc.**  
480 Neponset Street  
Bldg 12A  
Canton, MA 02021

Toll Free: 888-769-1246  
Phone: 781-828-0610  
Fax: 781-828-0542

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Web Site: [www.cellsciences.com](http://www.cellsciences.com)