

Mouse Anti-*Staphylococcus spp.* (Clone 17-5) mAb

Catalog No.	MON2008-1 MON2008-5	Quantity:	1.0 ml 5.00 ml
Description:	MON 2008 reacts with several <i>Staphylococcus</i> strains using ELISA. It does not react with several members of the Enterobacteriaceae.		
Concentration:	> 0.20 mg/ml		
Specificity:	<i>S. epidermis</i>		
Host:	Mouse		
Immunogen:	<i>S. epidermis</i>		
Isotype:	IgG3		
Clone:	17-5		
Formulation:	Tissue culture supernatant, RPMI, 0.02% sodium azide		
Applications:	ELISA, IF, agglutination		
Storage & Stability:	Store at 2-8°C for up to 1 year. Protect from light.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.