

Mouse Anti-Human Proteus spp (Clone 31-17) mAb

Catalog No.	MON2018-1 MON2018-5	Quantity:	1.0 ml 5.0 ml
Description:	The antibody reacts with bacteria in suspension and is therefore useful for detection of Proteus spp. in bacterial suspensions. Various systems can be used with this antibody to detect the presence of Proteus spp. MON 2017 reacts with approximately 50% of the Proteus spp. isolates tested.		
Concentration:	0.1 mg/ml		
Specificity:	Proteus spp.		
Host:	Mouse		
Isotype:	IgG2a		
Clone:	31-17		
Formulation:	Purified antibody in PBS containing 0.1% BSA and 0.02% sodium azide		
Storage & Stability:	Store at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.