

## FABP3

### Mouse Anti-Human H-FABP (Clone 67D3) mAb, Biotin

<b>Catalog No.</b>	MON2020B	<b>Quantity:</b>	0.5 ml
<b>Alternate Names:</b>	Fatty acid-binding protein, heart, Fatty acid-binding protein 3, Heart-type fatty acid-binding protein, H-FABP, Mammary-derived growth inhibitor, MDGI, Muscle fatty acid-binding protein, M-FABP		
<b>Description:</b>	<p>MON 2020 recognizes human heart-type fatty acid-binding protein (H-FABP) of both natural and recombinant origin. The H-FABP protein is derived from the human FABP3 gene. FABPs are small intracellular proteins (~13-14 kDa) with a high degree of tissue specificity that bind long chain fatty acids. They are abundantly present in various cell types and play an important role in the intracellular utilization of fatty acids, transport and metabolism. There are at least nine distinct types of FABP, each showing a specific pattern of tissue expression. Due to its small size, FABP leaks rapidly out of ischemically damaged necrotic cells leading to a rise in serum levels. Ischemically damaged tissues are characterized histologically by absence (or low presence) of FABP facilitating recognition of such areas. H-FABP is localized in the heart, skeletal and smooth muscle, mammary epithelial cells, aorta, distal tubules of the kidney, lung, brain, placenta, and ovary. It is also useful as marker for brain damage. Furthermore, this antibody is useful for the purification of H-FABP</p>		
<b>Concentration:</b>	0.1 mg/ml		
<b>Gene ID:</b>	P05413		
<b>Specificity:</b>	Binds to human heart Fatty Acid Binding protein. It shows no binding to human intestinal or liver FABP.		
<b>Immunogen:</b>	Purified human H-FABP		
<b>Isotype:</b>	Mouse IgG1		
<b>Clone:</b>	67D3		
<b>Conjugate:</b>	Biotin		
<b>Formulation:</b>	Biotinylated antibody solution in PBS, containing 0.02% sodium azide and 0.1% bovine serum albumin. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
<b>Applications:</b>	WB, IHC-F, IP		
<b>Storage &amp; Stability:</b>	Store in the dark at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

