

Cd8a, CD8b

Mouse Anti-Rat CD8 (Clone 15-11C5) mAb

Catalog No.	MON2023	Quantity:	1 ml
Alternate Names:	T-cell surface glycoprotein CD8 alpha and beta chains, CD8 antigen 32 kDa and 37 kDa chains, OX-8 membrane antigen, CD8a and CD8b		
Description:	MON 2023 recognizes rat CD8 antigen (MW 34 kDa and 39 kDa) and is reactive with all common inbred strains.		
Concentration:	0.1 mg/ml		
UniProt ID:	P07725, P05541		
Specificity:	Rat CD8		
Host:	Mouse		
Immunogen:	WAG/Rij spleen cells.		
Isotype:	IgG1		
Clone:	15-11C5		
Formulation:	Purified antibody solution in PBS, containing 0.1% BSA and 0.02% sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Applications:	IHC-F: suggested starting dilution 1:100 FC: suggested starting dilution 1:10		
Storage & Stability:	Store at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

