

Mouse Anti-Human Blood Group Antigen A1B (Clone HE-24) mAb

Catalog No.	MON2041	Quantity:	1.0 mL
Description:	MON 2041 distinguishes blood group A1B from A2B. The specificity of the antibody HE-24 was verified on >1000 samples of blood.		
Concentration:	1 mg/ml		
Specificity:	Human blood group antigen A1B		
Host:	Mouse		
Immunogen:	Glycoprotein fraction isolated from saliva of secretors with blood group A.		
Isotype:	IgM		
Clone:	HE-24		
Formulation:	Cell culture supernatant with 0.1% sodium azide		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

