

## Mouse Anti-Human Blood Group Antigen A (Clone H-193) mAb

<b>Catalog No.</b>	MON2044	<b>Quantity:</b>	1.0 mL
<b>Description:</b>	MON 2044 recognizes human blood group A (monofucosyl and difucosyl A antigens with chain types 1 and 2, A antigens with chain types 3, 4, 5, 6) and Forssman antigen.		
<b>Concentration:</b>	1 mg/ml		
<b>Specificity:</b>	Human blood group antigen A		
<b>Host:</b>	Mouse		
<b>Immunogen:</b>	Mixture of erythrocytes of group A1 and glycoprotein fraction isolated from saliva of secretors with blood group A.		
<b>Isotype:</b>	IgM		
<b>Clone:</b>	HE-193		
<b>Formulation:</b>	Cell culture supernatant with 0.1% sodium azide		
<b>Applications:</b>	FC, IHC		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

