

Mouse Anti-Human Blood Group Antigen A (Clone H-195) mAb

Catalog No.	MON2045	Quantity:	1.0 mL
Description:	MON 2045 recognizes human blood group A including weak variants A3, AX, A3B, AXB. The specificity of antibody HE-195 was confirmed by comparison of specificity and reactivity to standard reagent using >5.000 samples of blood.		
Concentration:	1 mg/ml		
Specificity:	Human blood group antigen A		
Host:	Mouse		
Immunogen:	Mixture of erythrocytes of group A1 and glycoprotein fraction isolated from saliva of secretors with blood group A.		
Isotype:	IgM		
Clone:	HE-195		
Formulation:	Cell culture supernatant with 0.1% sodium azide		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.