

Mouse Anti-Phosphotyrosine (Clone P-TYR-01) mAb

Catalog No.	MON2055	Quantity:	100 µg
Description:	MON 2055 detects tyrosine phosphorylation in activated cells and reacts with phosphotyrosine.		
Concentration:	1.0 mg/ml		
Specificity:	Phosphotyrosine		
Host:	Mouse		
Immunogen:	Phosphotyrosine conjugated to BSA		
Isotype:	IgG1		
Clone:	P-TYR-01		
Formulation:	Purified antibody in PBS, containing 0.1% bovine serum albumin and 15 mM sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Applications:	WB, ICC		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.