

## Mouse Anti-Human Hepatocyte (Clone EP-2) mAb

**Catalog No.** MON2059 **Quantity:** 100 µg

**Description:** MON 2059 stains human hepatocytes and the majority of human hepatocarcinomas. In frozen sections, it stains the cytoplasm of hepatic cells and may be used as a marker of the liver. Cell preparations of hepatocellular carcinoma biopsies or cell lines are found to bind EP-2 on the cell surface. Monoclonal antibodies to liver cell processes are useful in the identification of hepatic carcinomas and normal organ structures.

**Specificity:** Hepatocytes

**Host:** Mouse

**Isotype:** IgG2b

**Clone:** EP-2

**Formulation:** Purified antibody in PBS with 0.02% sodium azide. **Precaution:** Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.

**Applications:** IHC-F

**Storage & Stability:** Store at 2-8°C until expiration date.

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**