

## HLA-G

### Mouse Anti-Human HLA-G (Clone MEM-G/9) mAb, R-PE

<b>Catalog No.</b>	MON2075R	<b>Quantity:</b>	100 tests
<b>Description:</b>	MON 2075R reacts with native form of human HLA-G1 on the cell surface as well as with soluble HLAG5 isoform in its beta2 microglobulin associated form. Reactivity with HLA-G3 was also reported. Specific reaction was demonstrated with cells transfected with either full-length HLA-G and HLA-G1 complementary DNA.		
<b>Concentration:</b>	1 mg/ml		
<b>UniProt ID:</b>	P17693		
<b>Specificity:</b>	Human HLA-G		
<b>Immunogen:</b>	Recombinant human HLA-G refolded with beta 2-microglobulin and peptide		
<b>Host:</b>	Mouse		
<b>Isotype:</b>	IgG1		
<b>Clone:</b>	MEM-G/9		
<b>Conjugate:</b>	R-phycoerythrin		
<b>Formulation:</b>	Purified antibody in PBS containing 15 mM sodium azide		
<b>Applications:</b>	FC		
<b>Storage &amp; Stability:</b>	Store in the dark at 2-8°C until expiration date. <b>DO NOT FREEZE.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

