

SYK

Mouse Anti-Human Tyrosine Kinase SYK (Clone SYK-01) mAb

Catalog No.	MON2084	Quantity:	100 µg
Alternate Names:	Tyrosine-protein kinase SYK, Spleen tyrosine kinase, p72-Syk		
Description:	MON 2084 reacts with protein tyrosine kinase p72-Syk which is required for the transduction of signals through the B cell antigen receptor (BCR) and the high affinity IgE receptor (FcεRI).		
Concentration:	1 mg/ml		
UniProt ID:	P43405		
Purity:	>95% (by SDS-PAGE)		
Specificity:	Human Syk		
Host:	Mouse		
Immunogen:	Recombinant fragment (aa 5-360) of human Syk		
Immunoglobulin type:	IgG1		
Clone:	SYK-01		
Formulation:	Purified antibody in PBS containing 15 mM sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Cross-Reactivity:	Mouse, Rat		
Applications:	IHC, IP, WB, ICC		
Storage & Stability:	Store at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

