

HLA-G

Mouse Anti-Human HLA-G (Clone MEM-G/11) mAb, FITC

Catalog No.	MON2087F	Quantity:	100 tests
Alternate Names:	HLA class I histocompatibility antigen, alpha chain G, HLA G antigen, MHC class I antigen G		
Description:	MON 2087F recognizes HLA-G1 antigen. HLA-G belongs to the MHC Class I molecules (MHC Class Ib; nonclassical) and it is expressed on the surface of trophoblast cells.		
Concentration:	1 mg/ml		
UniProt ID:	P17693		
Specificity:	Human HLA-G		
Host:	Mouse		
Immunogen:	Recombinant human HLA-G refolded with beta 2-microglobulin and peptide		
Isotype:	IgG1		
Clone:	MEM-G/11		
Conjugate:	Fluorescein isothiocyanate (FITC)		
Formulation:	Purified antibody in PBS containing 15mM sodium azide		
Applications:	FC		
Storage & Stability:	Store at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

