

ALPP

Mouse Anti-Human Placental Alkaline Phosphatase (Clone H7E8) mAb

Catalog No.	MON2088	Quantity:	100 µg
Alternate Names:	Alkaline phosphatase, placental type, Alkaline phosphatase Regan isozyme, Placental alkaline phosphatase 1, PLAP-1		
Description:	MON 2088 recognizes throphoblast-derived membrane-associated enzyme alkaline phosphatase.		
Concentration:	1 mg/ml		
UniProt ID:	P05187		
Specificity:	Human placental alkaline phosphatase		
Host:	Mouse		
Immunogen:	Purified human placental alkaline phosphatase		
Immunoglobulin type:	IgG2a		
Clone:	H7E8		
Formulation:	Purified antibody in PBS containing 15mM sodium azide		
Applications:	WB, IHC, ELISA		
Storage & Stability:	Store at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

