

VIM

Mouse Anti-Human Vimentin (Clone VI-RE/1) mAb

Catalog No.	MON2091	Quantity:	100 ug
Alternate Names:	Vimentin		
Description:	MON 2091 reacts with human vimentin, a 57 kDa intermediate filament protein expressed on a wide variety of mesenchymal and mesodermal cell types.		
Concentration:	1 mg/ml		
UniProt ID:	P08670		
Specificity:	Human vimentin		
Host:	Mouse		
Immunogen:	Recombinant human vimentin, full length		
Isotype:	IgG1		
Clone:	VI-RE/1		
Formulation:	Purified antibody in PBS containing 15mM sodium azide		
Applications:	ELISA, FC, ICC, WB		
Storage & Stability:	Store at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

