

## KRT18

### Mouse Anti-Human Keratin-18 (Clone M9) mAb

<b>Catalog No.</b>	MON3006-1 MON3006-5	<b>Quantity:</b>	1 ml 5 ml
<b>Alternate Names:</b>	Keratin, type I cytoskeletal 18, Cell proliferation-inducing gene 46 protein, Cytokeratin-18, CK-18, Keratin-18, K18		
<b>Description:</b>	MON 3006 recognizes only keratin no. 18 (45 kD) and lacks cross reactivity with other cytokeratins. and no cross reactivity with GFAP, Keratin, neurofilament or Desmin detected.		
<b>Concentration:</b>	approximately 20 ug/ml		
<b>UniProt ID:</b>	P05783		
<b>Specificity:</b>	Human keratin 18 (45 kDa)		
<b>Host:</b>	Mouse		
<b>Immunogen:</b>	Cytokeratin isolated from human breast carcinoma cell line MCF-7		
<b>Isotype:</b>	IgG1		
<b>Clone:</b>	M9		
<b>Formulation:</b>	Cell culture supernatant in PBS containing 1% BSA, 0.1% sodium azide		
<b>Applications:</b>	<b>IHC-F:</b> suggested dilution 1:150 <b>IHC-P:</b> suggested dilution 1:10 <b>Western blot</b>		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

