

KRT8

Mouse Anti-Human Keratin 8 (Clone M20) mAb

Catalog No.	MON3016	Quantity:	1.0 mL
Alternate Names:	Keratin, type II cytoskeletal 8, Cytokeratin-8, CK-8, Keratin-8, K8, Type-II keratin Kb8		
Description:	MON 3016 stains almost all epithelial cells. In immunoblotting experiments the antibody recognizes only keratin no. 8 (52,2 kD) and does not cross react with other cytokeratins. Cross reactivity Vimentin Desmin, GFAP and Neurofilament not detectable.		
Concentration:	approximately 20 ug/ml		
UniProt ID:	P05787		
Specificity:	Human Keratin 8		
Host:	Mouse		
Immunogen:	Cytokeratin purified from the human breast carcinoma cell line MC-7		
Isotype:	IgG1		
Clone:	M20		
Formulation:	Cell culture supernatant in PBS containing 1% BSA, 0.1% sodium azide		
Applications:	IHC, WB		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

