

KRT13

Anti-Human Keratin 13 (Clone 1C7) mAb

Catalog No.	MON3017
Lot No:	1185001
Quantity:	1 ml
Alternate Names:	Keratin, type I cytoskeletal 13, Cytokeratin-13, CK-13, Keratin-13, K13
Description:	MON 3017 does not stain non-cornifying squamous epithelium. In immunoblotting experiments the antibody recognizes only keratin no. 13 (54kD) and does not show any cross reactivity with other cytokeratins. Cross reactivity Vimentin, Desmin, GFAP and Neurofilament not detectable.
UniProt ID:	P13646
Specificity:	Human Keratin 13
Host:	Mouse
Immunogen:	Cytokeratin purified from the human breast carcinoma cell line MC-7
Isotype:	IgG2a
Clone:	1C7
Formulation:	Cell culture supernatant in 0.15 M PBS containing 1% BSA and 0.1% sodium azide
Applications:	WB, IF, IP
Storage & Stability:	Store at 2-8°C up to 1 year.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

