

KRT13

Anti-Human Keratin 13 (Clone 2D7) mAb

Catalog No.	MON3018	Quantity:	1 ml
Alternate Names:	Keratin, type I cytoskeletal 13, Cytokeratin-13, CK-13, Keratin-13, K13		
Description:	MON 3018 does stain non-cornifying squamous epithelium. In immunoblotting experiments the antibody recognizes only keratin no. 13 (54kD) and does not show any cross reactivity with other keratins. Cross reactivity Vimentin, Desmin, GFAP and Neurofilament not detectable.		
Concentration:	approximately 20 ug/ml		
UniProt ID:	P13646		
Specificity:	Human Keratin 13		
Host:	Mouse		
Immunogen:	Cytokeratin purified from human esophagus		
Isotype:	IgG2b		
Clone:	2D7		
Formulation:	Cell culture supernatant in 0.15 M PBS containing 1% BSA and 0.1% sodium azide		
Applications:	WB, IHC		
Storage & Stability:	Store at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.