

## KRT18

### Mouse Anti-Human Keratin 18 (Clone DC-10) mAb

<b>Catalog No.</b>	MON3023	<b>Quantity:</b>	100 ug
<b>Alternate Names:</b>	Keratin, type I cytoskeletal 18, Cell proliferation-inducing gene 46 protein, Cytokeratin-18, CK-18, Keratin-18, K18		
<b>Description:</b>	MON 3023 reacts with the 45 kDa human keratin 18. The antigen isolated for immunization came from breast carcinoma cell line PMC-42. No cross reactivity with cow, pig, sheep, dog rabbit, rat or mouse.		
<b>Concentration:</b>	1 mg/ml		
<b>UniProt ID:</b>	P05783		
<b>Specificity:</b>	Human Keratin 18		
<b>Host:</b>	Mouse		
<b>Immunogen:</b>	Human breast carcinoma cell line PMC-42		
<b>Isotype:</b>	IgG1		
<b>Clone:</b>	DC-10		
<b>Formulation:</b>	Purified antibody in PBS containing 15 mM sodium azide		
<b>Applications:</b>	WB, IHC-P, FC, ICC, IP		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

