

KRT10

Mouse Anti-Human Keratin 10 (Clone DE-K10) mAb

Catalog No.	MON3025	Quantity:	1 mL
Alternate Names:	Keratin, type I cytoskeletal 10, Cytokeratin-10, CK-10, Keratin-10, K10		
Description:	MON 3025 specifically recognizes the 56.6 kD keratin K10 (formerly also designated cytokeratin 10) polypeptide in squamous cells and carcinomas from epidermis, lung, bladder, cervix, esophagus.		
UniProt ID:	P13645		
Specificity:	Human Keratin 10		
Host:	Mouse		
Immunogen:	Cytoskeletal preparation from human epidermis		
Isotype:	IgG1		
Clone:	DE-K10		
Formulation:	Cell culture supernatant containing 0.09 % sodium azide		
Applications:	Western blot: suggested dilution 1:100-1:250 IHC: suggested dilution 1:5-1:20		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

