

KRT5, KRT18

Mouse Anti-Human Keratin 5 and 18 (Clone C-50) mAb

Catalog No.	MON3029	Quantity:	100 ug
Alternate Names:	K5, K18, CK5, CK18, CYK5, CYK18, KRT5, KRT18		
Description:	MON 3029 reacts with Cytokeratin peptides 5 and 18 (58 and 45 kDa). Cytokeratins are a member of intermediate filaments subfamily represented in epithelial tissues.		
Concentration:	1 mg/ml		
Specificity:	Keratin 5 and 18		
Host:	Mouse		
Immunogen:	Cytoskeleton preparation from HeLa cells		
Isotype:	IgG1		
Clone:	C-50		
Formulation:	Purified antibody in PBS containing 15 mM sodium azide		
Cross-Reactivity:	Mouse, rat, hamster, pig, cow, sheep, and dog		
Applications:	IHC, IP, WB, ICC		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.