


KRT8

Mouse Anti-Human Keratin 8 (Clone C-51) mAb

Catalog No.	MON3032	Quantity:	100 ug
Alternate Names:	Keratin, type II cytoskeletal 8, Cytokeratin-8, CK-8, Keratin-8, K8, Type-II keratin Kb8		
Description:	MON 3032 reacts with Cytokeratin 8 (52.5 kDa). Cytokeratins are a member of intermediate filaments subfamily represented in epithelial tissues.		
Concentration:	1 mg/ml		
UniProt ID:	P05787		
Specificity:	Keratin 8		
Host:	Mouse		
Immunogen:	Cytoskeleton preparation from HeLa cells		
Isotype:	IgG1		
Clone:	C-51		
Formulation:	Purified antibody in PBS containing 15 mM sodium azide		
Cross-Reactivity:	Pig, bovine, sheep epithelial tissues		
Applications:	WB, ICC, IHC-P, IP		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

Western blot of human keratin 8



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.