

KRT19

Mouse Anti-Human Keratin 19 (Clone BA-17) mAb

Catalog No.	MON3033	Quantity:	100 ug
Alternate Names:	Keratin, type I cytoskeletal 19, Cytokeratin-19, CK-19, Keratin-19, K19		
Description:	MON 3033 reacts with cytokeratin 19 (40 kDa) in human tissue. Cytokeratin 19 is not expressed in hepatocytes; it is often co-expressed with cytokeratin 7.		
Concentration:	1 mg/ml		
Uniprot ID:	P08727		
Specificity:	Keratin 19		
Host:	Mouse		
Immunogen:	Mammary organoids		
Isotype:	IgG1		
Clone:	BA-17		
Formulation:	Purified antibody in PBS containing 15 mM sodium azide		
Cross-Reactivity:	Other mammalian species		
Applications:	WB, IP, ICC, IHC-P, FC		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

